

Trampoline - Rebound therapy model

230177UM



Details

Our Rebound Therapy trampoline is a 77 series school model trampoline with an UltraMesh bed and is approved by ReboundTherapy.org for use as a rebound therapy trampoline.

As well as Rebound Therapy, this trampoline can also be used for basic school trampolining. The trampoline is supplied with coverall frame pads.

This trampoline is supplied with either fixed height or hydraulically assisted "Safelift" roller stands

Dimensions:

- Frame size 4.6m x 2.75m x 1.0m.
- Folded size 2.75m x 2.0m x 0.8m.
- Bed size 3.6m x 1.8m

Height of trampoline when on stands:

- fixed height stands 2,000mm
- Safelift stands (high position) 2,010mm, Safelift stands (low position) 1,780mm

Options:

- Hydraulically assisted "Safelift" roller stands
- Safety coaching sides. These are extension pieces to the trampoline framework running the length of the trampoline on both sides. These, along with extra width coverall pads provide an area where a carer (or second carer) can stand without affecting the performance of the trampoline. We would recommend safety sides are specified to provide the most flexibility for rebound therapy. The upgrade price to add safety sides includes the extra cost for the additional width coverall pads to cover the wider sides.PLEASE NOTE: if you specify coaching sides, you cannot specify hydraulically operated Safelift roller stands as



these items are incompatible.


Bed:

The Rebound Therapy "UltraMesh" bed is a two layer bed with an UltraMesh material on the top which provides good grip for users in socks - which was specially developed in conjunction with the Rebound Therapy Organisation, and an under-layer for additional safety of genuine Permatron. The UltraMesh is available in a choice of colours:

- Black
- Grey
- Red
- Dark Blue
- · Light Blue
- Green

Coveralls:

As standard we supply contrasting coverall pads to the UltraMesh. As standard the colours are:






Coveralls are also available a wide range of <u>non-standard colours</u>. Please specify with your order the colour you require if it is not standard

Bed border:

As standard we sew a 50mm wide nylon border in BRIGHT GREEN around the perimeter of the bed to provide a visual clue as to where the bed ends and the pads start. This helps visually impaired people to discern where to keep their feet so as to avoid getting feet or hands under the pads and near the springs.



Alternative contrasting colours of border are available but BRIGHT GREEN will be used if no other colour is specified.

Steps:

We recommend our foam trampoline steps are used to provide easy and safe access onto the trampoline. These are available at additional cost - please see the product to the right:



61/11		
SKU	Options	Available Colours
230177UMGYNSL	Grey UltraMesh	N/A
230177UMGYSSF	Grey UltraMesh	N/A
230177UMGYNSF	Grey UltraMesh	N/A
230177UMRDNSF	Red UltraMesh	N/A
230177UMRDNSL	Red UltraMesh	N/A
230177UMRDSSF	Red UltraMesh	N/A
230177UMDBNSL	Dark Blue UltraMesh	N/A
230177UMDBSSF	Dark Blue UltraMesh	N/A
230177UMDBNSF	Dark Blue UltraMesh	N/A
230177UMGNNSL	Green UltraMesh	N/A
230177UMGNSSF	Green UltraMesh	N/A
230177UMGNNSF	Green UltraMesh	N/A
230177UMLBSSF	Light Blue UltraMesh	N/A
230177UMLBNSF	Light Blue UltraMesh	N/A
230177UMLBNSL	Light Blue UltraMesh	N/A



230177UMBKNSL	Black UltraMesh	N/A
230177UMBKSSF	Black UltraMesh	N/A
230177UMBKNSF	Black UltraMesh	N/A